*M.blue*<sup>®</sup>

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*"Everyone has a cross to bear - some more than others. The most important thing is not to give up and to live your life to its fullest."*

Sarah, nurse and hydrocephalus patient

## Foreword

Hydrocephalus is a disease that is frequently connected to the widest range of other existing problems. Some say that hydrocephalus is not a disease but the symptom of another underlying disorder.

Successful treatment of hydrocephalus only became possible in the 1950s. In a dramatic race for the life of his son Casey who had hydrocephalus, Philadelphia technician John D. Holder took only a few weeks to develop a silicone valve. The implantation of this valve in March 1956 marked a great step forward in the clinical treatment of this disease.

Christoph Miethke GmbH & Co. KG has built on the insights of over 50 years of valve treatment, and with the development of gravitation technology pioneered the development of valves that consistently consider the physical framework conditions in the drainage of cerebrospinal fluid, thus ensuring an intraventricular pressure in all body positions that matches that of a healthy person as closely as possible. These days, the primary aim of hydrocephalus therapy is no longer to prevent death or severe impairment as the result of the condition, but to enable patients to lead a life that is as normal as possible.

Your *M.blue* serves as a gravitational valve. The purpose of this handbook is to provide you and your family with a small insight into the treatment of hydrocephalus with the *M.blue* implant.

# Table of contents

Basics	20
Hydrocephalus	20
Operating principle of <i>M.blue</i>	22
Materials of <i>M.blue</i>	26
Patient pass	26
Treatment complications	26
Conduct after the operation	26
About the company	28

# Basics

The human brain is surrounded by a special fluid (scientific term: cerebrospinal fluid - CSF). Inside the head are several cerebral ventricles (scientific term: ventricles) where CSF is produced. These ventricles are interconnected with channels and constitute a complex drainage system.

CSF has the task of protecting the brain from mechanical damage, regulating intracranial pressure, maintaining moisture levels of the brain tissue and transporting metabolic products.

Fresh CSF is produced in the ventricles each day which is absorbed by the venous blood system. A healthy balance between the production and absorption of CSF is thus very important.

# Hydrocephalus

Hydrocephalus means that the balance between production and absorption of CSF is impaired. If more CSF is being produced than can be drained, the cerebral ventricles expand and a hydrocephalus develops ("hydro": Greek for water; "cephale": Greek for head). The ventricles exert pressure on the surrounding brain tissue and can cause neurological damage, some of which irreversible.

## Symptoms

Symptoms depend on the level of impairment and may include: headache, nausea, vomiting, coordination disorders, sleepiness up to the loss of consciousness. For children under the age of two, head circumference can increase significantly, because the bones in their skulls have not yet fused.

The causes of hydrocephalus are very varied and occur at any stage in life.


## Treatment methods

Even though there have been consistent efforts to find alternative treatment options to the implantation of a valve, e.g. by drug-based treatment or recently by minimally invasive surgical procedures, to this day there is not really an alternative in many cases to the implantation of a drainage system, known as a shunt.

In a hydrocephalus, intracranial pressure has to be reduced to within normal limits. To this end, in most cases a shunt system is implanted, creating a link between the cerebral ventricles and another bodily cavity, generally the abdomen, in order to drain excess CSF from the ventricles.

There are three types of shunt implantation:

1. Ventriculo-atrial (from the cerebral ventricle to the right atrium of the heart)
2. Ventriculo-peritoneal (from the cerebral ventricle to the abdominal cavity)
3. Lumbo-peritoneal (from the spinal canal to the abdominal cavity)


## Operating principle of *M.blue*

A proper understanding of the operating principle of *M.blue* requires a brief look at the operating principles of M.Blue requires a brief look at the physics that are involved in this type of drainage that are involved in this type of drainage: Creating a link between the cerebral ventricle and the abdominal cavity using a tube system might allow the CSF to drain without control, particularly when the patient stands upright, thus transporting too much of the fluid from the ventricles to the abdominal cavity. This is called overdrainage and must be prevented by all means and must be prevented using what ever means necessary. However, is it is entirely normal for a mobile human being to change position frequently: to stand up, walk, sit down or lie down. All this results in continuous physical changes in the drainage system.


In order to achieve controlled drainage of the cerebrospinal fluid, a valve such as *M.blue* is implanted together with the tube system. This combination of a tube system and a valve is known as a shunt system. Ideally, a valve such as *M.blue* ensures that CSF is only drained when the pressure inside the cerebral ventricles becomes too high. To that end, the valves have preset or adjustable opening pressures, i.e. a valve opens once a pre-defined intraventricular pressure has been reached, thus preventing over- or underdrainage. However, the optimum opening pressure in a prone position is not the same as in an upright one. This is the result of the physical laws of gravity (gravitation).

### Physical background

The cerebral ventricle and abdominal cavity connected by the shunt system can be exemplary compared with two connected fluid vessels. When these vessels are level with each other (as in a lying human being), the liquid contained in them tends to equalise. A low opening pressure setting in the shunt system ensures sufficient control in the lying position.


However, when there is a difference in height between the connected fluid vessels (as in a person in upright body position), the fluid contained in the vessels tends to follow gravity and run downwards. The height of a patient has a significant impact on the suction effect being created in the catheter in the upright position. The valve previously set to a low pressure now also opens at a low opening pressure in the upright position, thus ensuring that excessive cerebrospinal fluid drains while the patient is upright.

A single opening pressure would be a constant compromise between the prone and the upright body position, and would constitute a serious restriction to the quality of life.

**This means that different opening pressures are required for the various body positions.**

Adjustable gravitational unit

Sapphire ball

Tantalum weight

Rotor

Differential pressure unit

Sapphire ball

Micro-coil spring

## Mode of action of *M.blue*

The *M.blue* is a position-dependent valve. It consists of an adjustable gravitational unit and a fixed differential pressure unit.

## Horizontal position

In the horizontal position, the gravitational unit is always open and does not present any resistance. The opening pressure of *M.blue* in the horizontal position is characterised by the adjustable differential pressure unit. The differential pressure unit operates the same in all body positions and usually has a low opening pressure set to ensure that in case of increasing intraventricular pressure the valve facilitates the timely drainage of cerebrospinal fluid, particularly in the prone position.

## Vertical position

The gravitational unit makes use of the principle of gravity. As soon as the patient stands up, a tantalum weight is activated which additionally increases the valve opening pressure due to its gravitational force. The more upright the upper body of the patient is, the greater the opening pressure of the whole *M.blue*. This allows the valve to respond with an automatic opening pressure adjustment depending on the body position. Furthermore, the opening pressure of the gravitational unit can also be adjusted to individual patient needs after implantation through the skin, without requiring new surgery.

When using *M.blue plus* (*M.blue* with *proGAV 2.0*), the opening pressure of the differential pressure unit can also be adjusted after implantation through the skin, if necessary.

## Adjustment mechanism

*M.blue* is an adjustable valve. The integrated rotor brake ensures that an unintentional adjustment by externally acting magnetic fields is excluded. The adjustable gravitational unit of *M.blue* is equipped with a feedback mechanism. If external pressure is exerted on the valve an acoustic signal (a click sound) is audible, or a resistance can be felt, as soon as the rotor brake has been released due to the valve housing construction. In other words the valve shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released the rotor is adjustment-proof again. Although the click caused by releasing the rotor brake is easily audible before implantation, this can be considerably reduced after implantation and filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patients themselves or through the use of a stethoscope.

**Warning: Pressing the valve should be reserved for the attending physician.**

**Warning: The shunt system may be fitted with a pumpable reservoir. As too frequent pumping can result in excessive drainage (overdrainage) and thus unfavourable pressure conditions, this process should exclusively be performed by a doctor**

## Materials of *M.blue*

*M.blue* consists of high-quality materials that have been proven and standardised for use in implants; the main component is titanium. The stable housing ensures reduction of outside influences that could affect the valve function (such as external pressure). This results in a high level of functional assurance and thus a long operating life.

The catheters are made from silicone and are not made with natural rubber latex.

## Patient pass

Each *M.blue* is accompanied by a patient pass. This is completed by the treating physician and thus contains all of the important information for follow-up examinations.

## Treatment complications

The treatment of hydrocephalus with a shunt system is not always free from complications. As with every surgical intervention, there is a risk of infection. Unfortunately, problems can develop occasionally that could be directly or indirectly linked to the implanted valve. Such complications include blockages of the shunt system or undesired excessive or insufficient drainage of cerebrospinal fluid (over- or underdrainage).

## Conduct after the operation

Patients supplied with a shunt system are normally not restricted in their daily lives. You can safely travel by plane, do sports or go to work with *M.blue*. Only for the period immediately after the surgery, you should discuss with your doctor to what extent you should avoid physical exertion and for how long that would be advisable in your specific case.

**The following applies for potential MRI follow-up examinations:**  
*M.blue* is MR Conditional up to 3 Tesla.

## About the company

Christoph Miethke GmbH & Co. KG, based in Potsdam, has from 1992 concentrated on the development, production and distribution of innovative neurosurgical implants for the treatment of hydrocephalus.

Our daily aim is through the innovative development and production of high-quality reliable products to improve the life of hydrocephalus patients. One of our foremost task hereby in cooperation with hospitals, physicians and patients is primarily this: to understand. Again and again, we change perspectives and meet with our colleagues, partners and patients on an equal footing.

For detailed information about our company, please visit our website: [www.miethke.com](http://www.miethke.com).

More about the *M.blue* is available on our *M.blue* website: [www.thenewblue.com](http://www.thenewblue.com).

More about the operating principle of our valves is available here: [www.miethke.com/produkte](http://www.miethke.com/produkte)

Or use our app:

Apple: <https://itunes.apple.com/de/app/miethke/id450290015?mt=8>

Android: <https://play.google.com/store/apps/details?id=com.miethke.graviton>

You are also welcome to contact our Medical Devices Consultants:

+49 331 620 83-0 // [info@miethke.com](mailto:info@miethke.com) // [www.miethke.com](http://www.miethke.com)

Roland Schulz

Michaela Funk-Neubarth

Josefine Kehl

Thammo Weise

August von Hardenberg

Andreas Bunge

Jan MÜgel


Thoralf Knitter

For up-to-date information, visit our social media channels:

 [facebook.com/MiethkeDeutschland](https://facebook.com/MiethkeDeutschland)  [youtube.com/Miethke\\_International](https://youtube.com/Miethke_International)

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany