

Vasco[®] Protect blue

NON-STERILE ANTIMICROBIAL EXAMINATION GLOVES | DATA SHEET



B. Braun Melsungen AG confirms that Vasco[®] Protect blue gloves comply with the following standards and directives:

EC CERTIFICATES AND APPLIED STANDARDS

Medical Device Class I according to Medical Device Directive (MDD) 93/42/EEC

EN 455 1-4, ISO 11193-1, ASTM D6319

Personal Protective Equipment Category III according to Personal Protective Equipment Regulation (PPER) EU 2016/425

EN 420, EN 374, ISO 16523, ISO 16604, ASTM F1671, ASTM D6978

QUALITY CERTIFICATES

ISO 9001, ISO 13485

PERSONAL PROTECTIVE EQUIPMENT

Information and Declaration of Conformity according to PPER (EU) 2016/425:



www.bbraun.com/gloves-declarations-of-conformity

B. Braun Melsungen AG

A handwritten signature in blue ink, appearing to read "H. Gaudin".

Dr. Hans-Ulrich Gaudin
Head of Global Regulatory Affairs OPM Germany

Vasco® Protect blue

NON-STERILE ANTIMICROBIAL EXAMINATION GLOVES | REGULATORY INFORMATION

MEDICAL DEVICE INFORMATION

MDD 93/42/EEC (CLASS I), EN 455



PERSONAL PROTECTIVE EQUIPMENT INFORMATION

Tested in accordance with:

ISO 374-1/Type B



KPT

CE 2777 PPE Regulation (EU) 2016/425 (Cat. III); EN 420:2003+A1:2009

Code letter	Test chemical	EN 374-1:2016 Permeation level	EN 374-4:2013 Mean degradation
K	Sodium hydroxide 40%	Level 6	-2,8%
P	Hydrogen peroxide 30%	Level 2	42,3%
T	Formaldehyde 37%	Level 4	25,5%

Tested acc. to EN 16523-1:2015

Performance levels acc. EN 374-1:2016 +A1:2018	1	2	3	4	5	6
Measured breakthrough times (mins)	> 10	> 30	> 60	> 120	> 240	> 480

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. NOTE: Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation.

AQL 1.0

Resistance to bacteria and fungi	pass
Resistance to virus	pass

ISO 374-5:2016



VIRUS

EN 421:2010



Protection against particulate radioactive contamination.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

Vasco[®] Protect blue

NON-STERILE ANTIMICROBIAL EXAMINATION GLOVES | TECHNICAL DATA



SIZE	REF	GLOVE DIMENSIONS (EN 455)	
		Width of palm	Total length
	200/180 pcs.		
XS	9203315	≤ 80 mm	
S	9203327	80 ± 10 mm	
M	9203342	95 ± 10 mm	≥ 240 mm
L	9203357	110 ± 10 mm	
XL	9203368	≥ 110 mm	

PHYSICAL PROPERTIES

		Min. specification	Typical value
Wall thickness	Finger	0.050 mm	0.083 mm
	Palm	0.050 mm	0.060 mm
	Cuff	0.035 mm	0.051 mm
Force at break	During shelf life	6 N	7.31 N before ageing 8.11 N after ageing
Elongation at break	Before ageing	450%	558 %
	After ageing	400%	457 %
Tensile strength	Before ageing	18 MPa	44.10 MPa
	After ageing	16 MPa	48.12 MPa

GLOVE DESIGN

Colour	violet-blue
Shape	straight fingers, ambidextrous fitting
Cuff	rolled rim, regular cuff
Surface finish	finger textured
Inner glove surface	online chlorinated, powder-free


GLOVE MATERIAL

Nitrile butadiene rubber (NBR)	
Latex allergy risk	free of latex proteins

ACCELERATORS

Zn-dithiocarbamate	
Free of thiurames and mercaptobenzothiazoles MBT	

ANTIMICROBIAL PROPERTIES

 ANTIMICROBIAL GLOVES	Antimicrobial glove technology "AMG™" contributing to prevent bacterial cross-contamination
---	---

LOGISTIC INFORMATION

Dispenser pack dimensions	dispenser pack 200/180 pcs.	245 x 124 x 74 mm (L x W x H)
Transportation carton	10 dispenser packs	
Shelf life	3 years	
Storage conditions	store at room temperature, protect from dust, humidity, sun light and ozone	

Vasco[®] Protect blue

NON-STERILE ANTIMICROBIAL EXAMINATION GLOVES | BARRIER PROPERTIES – CHEMICALS



Tested by SATRA, UK in accordance with

EN 374-3: Protective gloves against chemicals and micro-organisms – Determination of resistance to permeation by chemicals.

CHEMICAL	CAS REGISTRY NO.	PERMEATION PERFORMANCE LEVEL	BREAKTHROUGH TIME
Chlorhexidine gluconate 4%	55-56-1	level 6	> 480 min
Ethanol 35%	64-17-5	level 1	> 10 min
Formaldehyde 35%	50-00-0	level 5	> 240 min
Formaldehyde 37%	50-00-0	level 4	> 120 min
Formalin 10%	50-00-0	level 6	> 480 min
Glutaraldehyde 1%	111-30-8	level 6	> 480 min
Glutaraldehyde 4%	111-30-8	level 6	> 480 min
Glycolic acid 2.5%	79-14-1	level 6	> 480 min
Hydrogen peroxide 3%	7722-84-1	level 5	> 240 min
Hydrogen peroxide 30%	7722-84-1	level 1	> 10 min
Nitric acid 65%	7697-37-2	not recommended	immediate
Sodium hydroxide 40%	1310-73-2	level 6	> 480 min
Sodium percarbonate 15%	15630-89-4	level 6	> 480 min

Vasco[®] Protect blue

NON-STERILE ANTIMICROBIAL EXAMINATION GLOVES | BARRIER PROPERTIES – CYTOSTATIC DRUGS



Tested by ARDL, USA in accordance with

ASTM D 6978: Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. Minimum detection rate < 0,01 µg/cm²/min

CLASSIFICATION

- Not suitable
- Suitable if changed before permeation breakthrough
- Suitable for prolonged use

CHEMOTHERAPY DRUG	MG/ML	CAS REGISTRY NO.	MIN BREAKTHROUGH DETECTION TIME
Carmustine	3.3	154-93-8	■ 16 min
Cisplatin	1.0	15663-27-1	■ > 240 min
Cyclophosphamide	20.0	6055-19-2	■ > 240 min
Dacarbazine	10.0	4342-03-4	■ > 240 min
Doxorubicin hydrochloride	2.0	25316-40-9	■ > 240 min
Etoposide	20.0	33419-42-0	■ > 240 min
Fluorouracil	50.0	51-21-8	■ > 240 min
Methotrexate	25.0	59-05-2	■ > 240 min
Mitomycin C	0.5	50-07-7	■ > 240 min
Paclitaxel (Taxol)	6.0	33069-62-4	■ > 240 min
Thio-Tepa	10.0	52-24-4	■ 9 min
Vincristine sulfate	1.0	2068-78-2	■ > 240 min