







## Infusionsvolumen | Volume à perfuser (ml)

	100		250		500		1000		2000		3000	
h	 ml	ml	 ml	ml	 ml	ml	 ml	ml	 ml	ml	 ml	ml
1	33	100	83	250								
2	17	50	42	125	83	250						
3	11	33	28	83	56	167	111	333				
4	8	25	21	63	42	125	83	250				
5	7	20	17	50	33	100	67	200	133	400		
6	6	17	14	42	28	83	56	167	111	333		
7	5	14	12	36	24	72	48	143	95	286	143	429
8	4	12	10	31	21	63	42	125	83	250	125	375
9		11	9	28	19	56	37	111	74	222	111	333
10	3	10	8	25	17	50	33	100	67	200	100	300
12		8	7	21	14	42	28	83	56	167	83	250
14		7	6	18	12	36	24	71	48	143	71	214
16	2	6	5	16	10	31	21	63	42	125	63	188
18				14	9	28	19	56	37	111	56	167
20		5	4	13	8	25	17	50	33	100	50	150
22				11		23	15	45	30	91	45	136
24	1	4	3	10	7	21	14	42	28	83	42	125

Infusionsdauer in Stunden = Durée de la perfusion en heures



= Tropfen/goutte(s)/min.

ml = ml/h



20 = 1 ml H<sub>2</sub>O



**B|BRAUN**  
SHARING EXPERTISE

Infusionstherapie  
von B. Braun

Thérapie par  
perfusion B. Braun

[bbraun.ch](http://bbraun.ch)

