



# Defix® | Technical Data Sheet

Single Use Kidney Bowls (non sterile)

		Defix®
REF		5385200
Dimensions	length	252 mm
	width	150 mm
	height	40 mm
	weight	13,4 g ± 5 %
Volume	gross volume	approx. 800 mL
	filling volume	approx. 500 mL
Material	type	polystyrene (PS), white
	thickness	0,40 mm ± 5 %
Logistics	packaging	box with 100 pieces
		carton with 4 boxes (400 pieces)
	shelf life	5 years
Sterilization		non sterile
Quality certificate		ISO 9001
Product classification		consumer product