

CIRCULAR STAPLER CURVED

ORDER INFORMATION

ARTICLE NO.	FIGURE	DESCRIPTION	COLOUR CODE	STAPLE QUANTITY	ANVIL OUTSIDE DIAMETER	LUMEN SIZE
CSC21		Conventional	Orange	16	21.0 mm	13.0 mm
CSC21-A		Air-Sealed	Orange	16	21.0 mm	13.0 mm
CSC25		Conventional	White	20	25.0 mm	17.0 mm
CSC25-A		Air-Sealed	White	20	25.0 mm	17.0 mm
CSC29		Conventional	Blue	24	29.0 mm	20.5 mm
CSC29-A		Air-Sealed	Blue	24	29.0 mm	20.5 mm
CSC33		Conventional	Green	28	33.0 mm	24.8 mm
CSC33-A		Air-Sealed	Green	28	33.0 mm	24.8 mm



GENERAL
OPEN SURGERY

CSC – CIRCULAR STAPLER CURVED

SINGLE USE INTERNAL STAPLER

Manufacturer acc. to directive 93/42/EEC:
Touchstone International Medical Science Co., Ltd. | 21 A Science Plaza
International Science Park | No. 1355 Jinjihu Avenue | Suzhou | China

Importeur:
B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Phone +34 93 586 6200 | Fax +34 93 588 1096 | www.bbraun.com

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

B. Braun Medical AG | Seesatz 17 | 6204 Sempach | Switzerland
Phone 0848 83 00 22 | Fax 0800 83 00 21 | sales-aesculap.bbmch@bbraun.com | www.bbraun.ch

The main product trademark "Aesculap" is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.




AS_05_0037_09.22

CIRCULAR STAPLER CURVED

SINGLE USE INTERNAL STAPLER

PRODUCT FEATURES


- RBS – Rotation Blade System. Easy & precise cutting.
- ASL – Built-in Automatic Safety Lock. Prevents accidental mis-firing.
- Bendos – Feeding tube re-entry pulley, protects anastomotic orifice.
- PSS – Petal Shaped Shaft. Easy insertion with firm attachment.
- Adjustable staple closure height ranging from 1 mm to 2.5 mm.



Before firing

After firing

DESCRIPTION



Bendos Feeding tube re-entry pulley

RBS Easy-cut Rotation Blade System

Slim trocar for efficient piercing