

ORDERING INFORMATION

HANDLE

Product Code	Shaft Length	Shaft Diameter
SELC-GS	60 mm	12 mm
SELC-GM	160 mm	12 mm
SELC-GL	260 mm	12 mm

LOADING UNIT

Product Code	Staple Line Length	Staple Height	Adapted Handle Product Code	Trocar Size	Color Code
SMG-30WR	30 mm	2.0/2.5/3.0 mm	SELC-GS SELC-GM SELC-GL	12 mm	White
SMG-30BR	30 mm	2.5/3.0/3.5 mm		12 mm	Blue
SMG-45WR	45 mm	2.0/2.5/3.0 mm		12 mm	White
SMG-45BR	45 mm	2.5/3.0/3.5 mm		12 mm	Blue
SMG-45DR	45 mm	3.0/3.5/3.8 mm		12 mm	Gold
SMG-45GR	45 mm	3.5/3.8/4.2 mm		12 mm	Green
SMT-45TR	45 mm	4.2/4.5/5.0 mm		15 mm	Black
SMG-60WR	60 mm	2.0/2.5/3.0 mm		12 mm	White
SMG-60BR	60 mm	2.5/3.0/3.5 mm		12 mm	Blue
SMG-60DR	60 mm	3.0/3.5/3.8 mm		12 mm	Gold
SMG-60GR	60 mm	3.5/3.8/4.2 mm		12 mm	Green
SMT-60TR	60 mm	4.2/4.5/5.0 mm		15 mm	Black

CARTRIDGE

Product Code	Staple Line Length	Staple Height	Adapted Loading Unit Product Code	Cartridge Color
SRG-30W	30	2.0/2.5/3.0 mm	SMG-30	White
SRG-30B	30	2.5/3.0/3.5 mm	SMG-30	Blue
SRG-45W	45	2.0/2.5/3.0 mm	SMG-45	White
SRG-45B	45	2.5/3.0/3.5 mm	SMG-45	Blue
SRG-45D	45	3.0/3.5/3.8 mm	SMG-45	Gold
SRG-45G	45	3.5/3.8/4.2 mm	SMG-45	Green
SRG-60W	60	2.0/2.5/3.0 mm	SMG-60	White
SRG-60B	60	2.5/3.0/3.5 mm	SMG-60	Blue
SRG-60D	60	3.0/3.5/3.8 mm	SMG-60	Gold
SRG-60G	60	3.5/3.8/4.2 mm	SMG-60	Green
SRT-45T	45	4.2/4.5/5.0 mm	SMT-45	Black
SRT-60T	60	4.2/4.5/5.0 mm	SMT-60	Black




SELC Bendos

ENDOSCOPIC LINEAR CUTTER AND
SINGLE USE LOADING UNIT

B. Braun Surgical, S.A.U. | Carretera de Terrassa, 121 | 08191 Rubi | Spain
Tel. +34 935 86 62 00 | +34 935 88 10 96 | www.bbraun.es

B. Braun Medical AG | Seesatz 17 | 6204 Sempach | Switzerland
Phone 0848 83 00 22 | Fax 0800 83 00 21 | sales-aesculap.bbmch@bbraun.com | www.bbraun.ch

Information according to Regulation (EU) 2017/745 and Council Directive 93/42/EEC respectively:

 **Manufacturer SELC:**
Touchstone International Medical Science Co., Ltd. | 278 Dongping Street | Suzhou Industrial Park | 215123 Suzhou | China
CE 0123

 **Importer of Product Y and Z:**
B. Braun Surgical, S.A.U. | Carretera de Terrassa, 121 | 08191 Rubi | Spain

 **Authorized Representative**
Shanghai International Holding Corp. GmbH (Europe) | Eiffestraße 80 | 20537 Hamburg | Germany

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form, partial or otherwise, is not permitted.

SELC Bendos

ENDOSCOPIC LINEAR CUTTER AND SINGLE USE LOADING UNIT



PRODUCT FEATURES AND BENEFITS



➤ **ONE** Platform
TWO Options of Reloads

➤ **FREEDOM** of Choice

Multiple options of reloads with interchangeable cartridges and different staple heights



2

Smart Combinations




3

TGC - Tissue Gap Control



Different Staple Heights



4