



LAPAROSCOPIC SURGERY

# AESCULAP® ADTEC® BIPOLAR


THE BLUE (R)EVOLUTION

# AESCULAP® ADTEC® BIPOLAR

THE BLUE (R) EVOLUTION







This product was developed through the use of materials – ceramic, steel and plastic – and combining these components using well established manufacturing technologies.

There is always room for improvement. So AESCULAP® has systematically refined these instruments.

Why not join the **AdTec**® **bipolar** and experience these exciting bipolar endoscopic instruments for yourself!









### **CONCEPT**

**AdTec® bipolar** stands for a bipolar instrument technology. Proofed application materials stand for stability, high tactility and, of course, high quality.

Patent: DE 196 08 716

US 5 853 412

### RELIABILITY

#### Reliability up to the tip

- The design enables good coagulation results.
- Reliable closure of even arterial vessels.
- Reduced risk of "slipping effect" when grasping and preparing.
- No current leakage thanks to bipolar technology.
- Reliable grasping, cutting and precise dissection of delicate structures.

### **TRANSMISSION**

The feeling of the touch is transmitted back to the user, letting him decide how firmly he wants to grasp the tissue.

### **UNIVERSAL**

#### Modular design

 Simple concept, comparable to the monopolar instruments: the same features as the monopolar instruments.

#### Universal use

- Can be used in laparoscopic indications in gynaecology, urology, general and vascular surgery.
- For tissue and delicate structures.

# AESCULAP® ADTEC® BIPOLAR

**COST EFFECTIVENESS** 


The **AdTec**® **bipolar** products are reusable and their modular design allows for inexpensive replacement of individual components as needed.











		Instrument complete	Jaw inserts	Insulated outer tube	Handle
METZENBAUM Scissors	Working length 220 mm	PM409R	PM439R	PM970R	PM450R
	Working length 310 mm	PM400R	PM430R	PM973R	PM450R
	Working length 420 mm	PM404R	PM434R	PM981R	PM450R
MINI METZENBAUM Scissors	Working length 220 mm	PM414R	PM444R	PM970R	PM450R
	Working length 310 mm	PM407R	PM437R	PM973R	PM450R
MARYLAND Fixation forceps	Working length 220 mm	PM410R	PM440R	PM970R	PM450R
	Working length 310 mm	PM401R	PM431R	PM973R	PM450R
	Working length 420 mm	PM405R	PM435R	PM981R	PM450R
MARYLAND Fixation forceps	Working length 220 mm	PM413R	PM443R	PM970R	PM450R
fenestrated	Working length 310 mm	PM408R	PM438R	PM973R	PM450R
	Working length 420 mm	PM415R	PM445R	PM981R	PM450R
Micro fixation forceps	Working length 310 mm	PM402R	PM432R	PM973R	PM450R
Macro fixation forceps	Working length 220 mm		PM442R		PM450R
, ,	Working length 310 mm	PM403R	PM433R	PM973R	PM450R
	Working length 420 mm	PM406R	PM436R	PM981R	PM450R
Protecting cap		PM433220		-	

# AESCULAP® ADTEC® BIPOLAR

REUSABLE AND DISPOSABLE CABLES

Reusable cables, 4 m length	Equipment side		Instrument side - un	Instrument side - universal connector	
GN132	ERBE int. Dolly EMC	22 mm			Aesculap round pin Aesculap tab US round pin
GN133	Aesculap BOWA Lamidey Söring Valleylab EMC	28.6 mm			Aesculap round pin Aesculap tab US round pin
GN134	Codman Synergy	38 mm	1		Aesculap round pin Aesculap tab US round pin
GN131	Erbe Storz	Ø 8 mm Ø 4 mm	(0)		Aesculap round pin Aesculap tab US round pin
GN130	Aesculap ACMI Berchtold Martin Wolf	Ø 5 mm Ø 2 mm			Aesculap round pin Aesculap tab US round pin
isposable cables, 3.7 length  GK290SU	ERBE int. Dolly EMC	22 mm			(II) (•)  Aesculap tab  US round pin
GK291SU	Aesculap BOWA Lamidey Söring Valleylab EMC	28.6 mm			Aesculap tab US round pin
GK292SU	Codman CMC I CMC II*	38 mm	INI		Aesculap tab US round pin
		+	GN088	Bipolar hand switch for AdTeo for use with HF generator Nel	•
>				For further cables and adap C30481.	ters see our brochui

Article No.	Description		
JG600	STERILIT		
	AESCULAP Oil Spray in 300 ml can (ozone-friendly), for maintenance of instruments prior to sterilization. Physiologically inert acc. to DAB and		
	§ 31 LMBG (German Pharmacopoeia).		
JG598	STERILIT		
	AESCULAP Instrument Oil for box locks, ratchets etc. prior to sterilization. Physiologically inert acc. to DAB and § 31 LMBG (German Pharmacopoeia).		





# ADTEC® SINGLE USE BIPOLAR SINGLE USE INSTRUMENTS

- Good cutting performance
- Jaw parts feature ceramic insulation
- Insulated to the tips

All Aesculap endoscopy products on our Internet homepage:



...just a click away!

### AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com