

OVER- AND UNDERFILLING OF EASYPUMP® II

SHORT TERM INFUSION THERAPY

PRODUCT SPECIFICATIONS

ST-Variants										
ST-Variants	100-0.5	250-0.5	50-1	100-1	250-1	250-1.5	400-2	500-2	100-2	400-4
Code No. (REF)*	4540040	4540042	4540044	4540046	4540048	4540050	4540052	4540054	4540056	4540058
Nominal Flow Rate (ml/h)	200	500	50	100	250	175	200	250	50	100
Nominal Fill Volume (ml)	100	250	50	100	250	250	400	500	100	400
Minimal Fill Volume (ml)	50	125	30	50	135	135	240	240	50	240
Maximal Fill Volume (ml)	125	295	65	125	335	335	560	560	125	560
Maximum Retained Volume (ml)	≤ 3	≤ 8	≤ 2	≤ 3	≤ 8	≤ 8	≤ 10	≤ 10	≤ 3	≤ 10

STORAGE

Time to reach room temperature (approx. +23 °C to ±2 °C)										
Hours from Refrigerator (approx. +2 °C to +8 °C)	6	12	6	6	12	12	12	12	6	12
Hours from Freezer (approx. -18 °C)	12	18	12	12	18	18	18	18	12	18

OVER- AND UNDERFILLING DATA

Approximate delivery time										
Hours (hh:mm)	Fill volume (ml)									
00:15 h	50	125	-	-	-	-	-	-	-	-
00:30 h	100	250	-	50	-	-	-	-	-	-
00:45 h	-	-	38	75	188	-	-	-	-	-
01:00 h	-	-	50	100	250	175	-	250	50	-
01:15 h	-	-	63	125	-	219	-	313	63	-
01:30 h	-	-	-	-	-	263	300	375	75	-
01:45 h	-	-	-	-	-	-	350	438	88	-
02:00 h	-	-	-	-	-	-	400	500	100	-
02:15 h	-	-	-	-	-	-	450	-	113	-
02:30 h	-	-	-	-	-	-	500	-	125	250
02:45 h	-	-	-	-	-	-	550	-	-	275
03:00 h	-	-	-	-	-	-	-	-	-	300
03:30 h	-	-	-	-	-	-	-	-	-	350
04:00 h	-	-	-	-	-	-	-	-	-	400
04:30 h	-	-	-	-	-	-	-	-	-	450
05:00 h	-	-	-	-	-	-	-	-	-	500
05:30 h	-	-	-	-	-	-	-	-	-	550

Note: In order to deliver the Nominal Fill Volume, please add the Maximum Retained Volume:

Example LT 60-12: 60 ml + 2 ml = 62 ml (Nominal Fill Volume + Maximum Retained Volume = Total Filling Volume)

OVER- AND UNDERFILLING OF EASYPUMP® II LONG TERM INFUSION THERAPY

PRODUCT SPECIFICATIONS

LT-Variants														
LT-Variants	60-12	500-12.5	80-16	125-25	270-27	60-30	120-30	400-40	100-50	270-54	400-80	270-68	400-100	270-135
Code No. (REF)	4540002	4540003	4540004	4540006	4540008	4540010	4540012	4540014	4540016	4540018	4540022	4540026	4540028	4540032
Nominal Flow Rate (ml/h)	5	40	5	5	10	2	4	10	2	5	5	4	4	2
Nominal Fill Volume (ml)	60	500	80	125	270	60	120	400	100	270	400	270	400	270
Minimal Fill Volume (ml)	30	240	50	60	120	30	60	240	50	135	240	135	240	135
Maximal Fill Volume (ml)	65	560	125	125	335	65	125	560	125	335	560	335	560	335
Maximum Retained Volume (ml)	≤2	≤10	≤3	≤3	≤8	≤2	≤3	≤10	≤3	≤8	≤10	≤8	≤10	≤8

STORAGE

Time to reach room temperature (approx. +23 °C to ±2 °C)														
Hours from Refrigerator (approx. +2 °C to +8 °C)	6	12	6	12	12	6	12	12	6	12	12	12	12	12
Hours from Freezer (approx. -18 °C)	12	18	12	18	18	12	18	18	12	18	18	18	18	18

OVER- AND UNDERFILLING DATA

Approximate delivery time														
Hours (hh:mm)	Days	Fill volume (ml)												
6 h	-	-	240	-	-	-	-	-	-	-	-	-	-	-
8 h	-	40	320	-	-	-	-	-	-	-	-	-	-	-
10 h	-	50	400	50	-	-	-	-	-	-	-	-	-	-
12 h	-	60	480	60	60	120	-	-	-	-	-	-	-	-
18 h	-	-	-	90	90	180	36	72	-	-	-	-	-	-
24 h	1 d	-	-	120	120	240	48	96	240	-	-	-	-	-
30 h	-	-	-	-	-	-	60	120	300	60	150	-	-	-
48 h	2 d	-	-	-	-	-	-	-	480	96	240	240	-	-
60 h	2.5 d	-	-	-	-	-	-	-	-	120	-	300	240	240
72 h	3 d	-	-	-	-	-	-	-	-	-	-	360	288	288
96 h	4 d	-	-	-	-	-	-	-	-	-	-	480	-	384
120 h	5 d	-	-	-	-	-	-	-	-	-	-	-	-	480
144 h	6 d	-	-	-	-	-	-	-	-	-	-	-	-	-
168 h	7 d	-	-	-	-	-	-	-	-	-	-	-	-	336
192 h	8 d	-	-	-	-	-	-	-	-	-	-	-	-	-
216 h	9 d	-	-	-	-	-	-	-	-	-	-	-	-	-
240 h	10 d	-	-	-	-	-	-	-	-	-	-	-	-	-
264 h	11 d	-	-	-	-	-	-	-	-	-	-	-	-	-
288 h	12 d	-	-	-	-	-	-	-	-	-	-	-	-	-

Note: In order to deliver the Nominal Fill Volume, please add the Maximum Retained Volume:

Example LT 60 - 12: 60 ml + 2 ml = 62 ml (Nominal Fill Volume + Maximum Retained Volume = Total Filling Volume)

B. Braun Medical AG | Seesatz 17 | 6204 Sempach | www.bbraun.ch
 Hospital Market: Tel. 0848 83 00 44 | Fax 0800 83 00 43 | sales-hospital.bbmch@bbraun.com
 Out Patient Market: Tel. 0848 83 00 33 | Fax 0800 83 00 32 | sales-opm.bbmch@bbraun.com